

www.biotium.com Revised: February 7, 2019

# **Product Information**

## Pluronic® F-127, 20% in DMSO

Catalog Number: 59004

Unit Size: 1 mL

**Properties** 

Molecular Weight: ~ 12500 **CAS Number:** 9003-11-6 Color & Form: White liquid

#### Storage and Handling

Store the solution at room temperature. DO NOT refrigerate or freeze. Product is stable for at least 12 months from date of receipt when stored as recommended.

Pluronic® F-127, 20% in DMSO may solidify into a gel during storage, especially at low temperature. This is normal and does not affect the product, but the solution must be in liquid form before use. Heat the vial to 50-65°C for 5-10 minutes and vortex periodically until it has formed a clear liquid.

#### **Product Description**

Pluronic® F-127 is a nonionic detergent useful for solubilizing hydrophobic molecules in aqueous solutions. Pluronic® F-127 is commonly used to solubilize hydrophobic AM ester forms of fluorescent indicator dyes for calcium and other ions for loading of the dyes into cells (see References). It also can be used to solubilize other hydrophobic dyes or compounds to facilitate cell penetration.

### **Assay Protocol**

The following is a typical procedure for loading AM esters of ion indicators into cells. Optimal conditions for cell loading may vary for different cell types or compounds.

- 1. Dissolve the dye or AM ester in anhydrous DMSO at 1-5 mM (or ~1000X to 500X the final desired dye concentration).
- 2. Immediately before use, mix equal volumes of the dye or AM ester stock solution and Pluronic® F-127, 20% in DMSO in a microcentrifuge tube.
- 3. Add the solution from step 2 to cell culture medium or buffer to achieve a final AM ester concentration of 1 uM to 10 uM.
- 4. Remove culture medium from the cells and add the medium or buffer containing diluted dye or AM ester with Pluronic® F-127.
- 5. Incubate cells at room temperature or 37°C for 10 minutes to 1 hour or longer.

#### References

- 1) J. Membrane Biol. 19, 1 (1974).
- 2) Science 233, 886 (1986).
- 3) J. Biol Chem 262, 12801 (1987).
- 4) J. Biol. Chem. 265, 19543 (1990).
- 5) Methods Enzymol. 302, 341 (1999).
- 6) Methods Enzymol. 307, 441(1999).

#### Related Products

Cat. No.	Product
90082	DMSO, Anhydrous
59000	Pluronic® F-127
59005	Pluronic® F-127, 10% in H <sub>2</sub> O
59100	Calcium Calibration Buffer Kit
59001	A-23187, Free Acid
59006	4-Bromo A-23187, Free Acid
59007	Ionomycin, Calcium Salt
59002	EDC (EDAC)
41024-4L	Water, Ultrapure Molecular Biology Grade
22023	Paraformaldehyde, 4% in PBS, Ready-to-Use Fixative
50015	Fluo-3, AM Ester, 1 mM in Anhydrous DMSO
50018	Fluo-4, AM Ester
50034	Fura-2, AM Ester
50038	Furaptra, AM Ester
50044	Indo-1, AM Ester
50023	Rhod-2, AM Ester
50025	Rhod-590, AM Ester
50000	BAPTA, AM Ester
50005	5,5'-Difluoro BAPTA, AM Ester
50007	5',5'-Dimethyl BAPTA, AM Ester

Please visit our website at www.biotium.com for information on our life science research products, including live cell stains for nuclei, mitochondria, and other organelles, apoptosis probes, fluorescent ion indicators, and labeled antibodies and other conjugates for cell biology research.

Pluronic is a registered trademark of BASF.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.